

Configuration

Canto II

Laser	Detector	Filter	Fluorochromes	Proteins
Blue (488nm)	Green	530/30	FITC, AF488, CFSE	EGFP, EYFP, mCitrine, Venus
	Orange	585/42	PE,	DsRed
	Red	695/40	PerCP, PerCP-Cy5.5, 7AAD, PI	
	Infrared	780/60	PE-Cy7, PE-Vio780	

Red (633nm)	Red	660/20	APC, Cy5, AF647	mPlum
	Infrared	780/60	APC-Cy7, APC-H7, APC-Vio780	

Violet (405nm)	Blue	450/50	PacificBlue, AF405, BV421, eF450, Horizon V450, VioBlue	Cerulean, ECFP
	Green	510/50	PacificOrange, Horizon V500, BV510, AmCyan, Q-Dot525, VioGreen	